

SAF-T-LOG[®]

Paperless HACCP in 15 Minutes



ThermoWorks

Simplified HACCP

Set high and low targets for each checklist item.
Corrective actions are displayed automatically if limits are breached.

Interchangeable type K probe.

Real-time Pass/Fail results.

Checklists in °F or °C.

Backlight for low-light conditions.


Store readings with date/time, operator, Pass/Fail, temp and corrective action.

Displays temperature checklist item or Yes/No question.

Sealed waterproof enclosure.

Quickly scroll through food and Yes/No items.

Super-Fast® reduced-tip probe reads temp in 3-4 seconds.





Scan to see the Saf-T-Log virtual product demo.

“It literally took 15 minutes to create a checklist and start using the Saf-T-Log.” — Quality Manager, San Gennaro Foods

Easy Setup

Step 1 → **CREATE** **Step 2** → **RECORD** **Step 3** → **SAVE**



Food Town Grocery
123 Park Drive
Naperville, IL 60540
Phone 312-555-5555
Fax 312-777-7777

Captured: 24-02-2012 14:20:30
Configuration: Morning Deli
SAF-T name: NF-214 (Supervisor)
Location: Deli Counter - Northside

ITEM	USER	DATE/TIME	LIMITS (°F)	READING (°F)	NOTE
Chicken Strips	Mike W.	24-02-2012 14:20:06	Hi: 150.0 Lo: 138.0	141.8	
Roasted Chicken	Mike W.	24-02-2012 14:20:08	Hi: 150.0 Lo: 138.0	135.6	Reheated product
Orange Chicken	Mike W.	24-02-2012 14:20:10	Hi: 165.0 Lo: 140.0	134.3	Reheated to 165°F
Potato Wedges	Mike W.	24-02-2012 14:20:11	Hi: 150.0 Lo: 138.0	140.2	
Nuggets	Rachel S.	24-02-2012 14:20:14	Hi: 150.0 Lo: 138.0	146.3	

Load the Saf-T-Log software* onto your computer and start building a custom checklist. Connect Saf-T-Log (via USB) to the computer and download the list you just created.

Scroll to the food item, insert the probe and press record. Scroll to the next item or question and repeat. It's really that easy!

Reconnect Saf-T-Log to your computer and temperature readings automatically download. Store them electronically. You can also email them to a centralized location. You may even choose to print - if you hate trees.

Incredible Value

Features

- List temperature and YES/NO questions
- Select corrective actions
- See real-time PASS/FAIL results
- Store date, time and user with each reading
- Automate data archiving and report creation
- Add as many as 300 items per checklist
- Download up to seven checklists to a single Saf-T-Log with an unlimited number available on your PC

Specifications

Range	-148 to 2500°F (-100 to 1372°C)
Resolution	0.1° to 600°F (300°C) thereafter 1°
Accuracy	±0.7°F (±0.4°C)
Probe	Type K Thermocouple
Water Resistance	IP66/67 rated
Connection Type	USB 2.0
Calibration	NIST-Traceable Certificate included

Quick & Simple HACCP

Saf-T-Log is a total innovation from ThermoWorks. No other handheld system records data this intuitively at anywhere near this price. In less than 15 minutes, you can set up your custom daily checklists, download them to the Saf-T-Log and start recording temperatures. In less than 15 minutes, you can effectively eliminate the need for paper log books, clip boards and tattered binders.

Streamline your HACCP recording, consolidate your daily checklists into one easy-to-use device and make temperature monitoring fun.

Ask an expert! Call us with any questions about the Saf-T-Log at: 1.800.393.6434, or email us at: techsupport@thermoworks.com

** There are NO annual fees or start-up costs associated with the Saf-T-Log. All software is included. Saf-T-Log includes a 1-year warranty against manufacturing defects.*

“This durable, accurate thermometer is easy to use, helps your bottom line and saves trees by eliminating paper logs. It makes temperature monitoring fun and eliminates hated paperwork.”

— **Senior Director, Global Food Safety - Marriott**

Cloud Accessible Reports

Saf-T-Log makes it easy to get the reports you need regardless of where you are. Using third-party cloud storage solutions like DropBox, you can access your Saf-T-Log reports anywhere. Now it's easier than ever for regional managers or corporate leaders to view completed checklists - immediately!

Accessories



830-257

Protective Rubber Boot



830-259

Protective Boot with Probe Holder

Go Paperless!

Only \$199

To Order

292-701	Saf-T-Log Recording Thermocouple Meter	\$199
Y-292-701	Saf-T-Log Kit w/313-158 fast resp. probe (\$25 savings)	\$224
313-158	Fast Response Probe	\$51
830-257	Protective Rubber Boot	\$14
830-259	Protective Boot with Probe Holder	\$14

To Order:

Visit: www.thermoworks.com/saf-t-log

or call: 1-800-393-6434

Saf-T-Log[®] Paperless HACCP in 15 Minutes




ThermoWorks

1762 W. 20 S. #100
Lindon, UT 84042



“Not only is the thermometer convenient, it saves a lot of time when conducting a line check. Our staff is now excited to conduct them as opposed to before when they moaned how long it took. Now it's a fun thing to do!”

— Regional Manager • Moe's Southwest Grill



Order the Saf-T-Log Kit
and save \$25

Order Today...



Go Paperless Tomorrow!