



ThermoWorks

WAND Blue™

Non-Contact, Digital Forehead Thermometer
with *Bluetooth*® Wireless Technology



No-Touch Infrared Technology
For safety and hygiene.

FDA-Cleared
For use in adults and children.

Forehead Use
Displays "oral equivalent"
temperature.

Integrate
Use our app or integrate into
your existing apps.

Quick 1 Second Readings
Instant results - no waiting.
Beeps when done.

Simple to Use
Point the thermometer at the forehead
and press Start. It's that easy!

Pairs Quickly
Pairs quickly and easily with
Bluetooth to your smart device.

High Temperature Alert
Audible alert when temperatures
are $\geq 99.5^{\circ}\text{F}$ (37.5°C)

Mute Setting
Disable audible alert that could
disturb a sleeping child.

Memory
Stores the 25 most recent readings
for easy recall or see in app.

Bright LED backlight
Makes temperatures easy to read,
even in the dark.



**Non-contact, 1-second readings for
quick, safe screening.**



**FDA-cleared for safety & accuracy for kids and adults.
ThermoWorks Wand Thermometer measures infrared energy
radiated from the skin at the center of the forehead area. This
captured energy is collected through the lens and converted to a
body temperature value (displays oral equivalent temperature).*

*It's helpful to know each individual's normal temperature when
they are well. This is the best way to accurately diagnose a fever.
To get a reference value, we suggest taking multiple readings (5
to 10) over one minute, taking the highest number to determine
"normal temperature".*

The Bluetooth® word mark and logos are registered trademarks owned
by Bluetooth SIG, Inc. and any use of such marks by ThermoWorks, Inc. is
under license.